





Parts Required: Padded Thigh Belt (color and size specific)

Tools Used: No tools required

Step	Execution	Figures
1	<p>Locate the front cross piece bar of the seat frame. Starting from under the seat frame, pass one end of the thigh belt between the second and third strap from the front of the seat.</p> <p><i>Note: The padded thigh belt does not accommodate the full vinyl seat cover, as there is no "V" notch in the seat pad.</i></p>	
	<p>For Pedal Rocker attachment, locate the front frame piece attached to the seat frame (resembles a handle)</p>	
2	<p>Pass the other end of the Thigh Belt through the same spot in the chair seat. Ensure that the center strap of the Thigh Belt is looped around the front cross piece of the seat frame.</p>	
3	<p>One padded section of the Thigh Belt passes over each thigh. Connect the buckles behind the chair. Do not fasten around the cylinder or cable behind the chair. Tighten the straps as necessary.</p>	

Auto Style buckle also available. Repeat steps 1-4



For more information, please contact Broda Customer Care at 1-800-668-0637