

Additional Positioning Padding Package (APP)



*Shown on Synthesis
Positioning Wheelchair*

RECOMMENDED FOR:

Use for postural alignment with Broda positioning wheelchairs to give patients extra comfort.

PACKAGE INCLUDES:


- Head support with elevated edges
- Independently adjustable pillow
- Contoured back with lumbar support for postural alignment and centering
- High-specification contoured seat cushion for pelvic immersion
- Calf support for centering of lower extremities
- Foot support with elevated edges



Easy to install and attaches securely to chair



High performance foam covered with removable waterproof cover



Works in harmony with Comfort Tension Seating[®]



Enhanced body centering, safety, and pressure redistribution

Slim Posture Padding (SPP) Package

RECOMMENDED FOR:

Use for postural alignment with Broda pedal wheelchairs to give patients extra comfort with accommodations for easy foot propulsion.

PACKAGE INCLUDES:

- Back extender cushion for tall back chairs
- Back cushion with lumbar and lateral torso support
- Contoured seat cushion for offloading pressure



Shown on the Encore Pedal Wheelchair




Low profile and unobtrusive cushion



Lumbar support for spinal alignment and improved posture



Works in harmony with Comfort Tension Seating[®]



Gel-infused visco swirl reduces heat buildup and provides comfort.