

BRODA[®]

Product Lookbook

Essential Comfort



BRODA[®]

brodaseating.com
1-844-552-7632

TABLE OF CONTENTS

3	<i>TYPES OF CHAIRS</i>
6	<i>SEATING & POSITIONING ESSENTIALS</i>
8	<i>POSITIONING</i>
13	<i>PEDAL</i>
16	<i>REHABILITATION (E1161)</i>
21	<i>TRANSPORT</i>
23	<i>HYGIENE</i>
27	<i>BEHAVIORAL</i>
29	<i>PACKAGES & ACCESSORIES</i>
32	<i>CHAIR SPECIFICATIONS</i>
35	<i>UPHOLSTERY & FRAME COLORS</i>
39	<i>CONFIGURE YOUR WHEELCHAIR</i>

A photograph showing a woman in blue scrubs adjusting a wheelchair for a woman in a yellow shirt sitting in it. They are outdoors under a tree. The wheelchair has 'BRODA' written on the seat. A yellow banner is overlaid on the left side of the image.

*Essential
Comfort*

PRODUCT CATEGORIES

For over four decades, Broda has been bringing comfortable seating solutions to patients and caregivers in long-term and complex care settings. Our full line of tilt & recline positioning wheelchairs, pedal wheelchairs, rehab wheelchairs, shower commodes, and gliders provides essential comfort and increased mobility for wheelchair users all over the world.

POSITIONING



Tilt and recline functionality provide the highest level of comfort and care

PEDAL



Broda wheelchairs that enhance an individual's ability to safely move and be independent

REHABILITATION (E1161)



Durable, versatile wheelchairs give support in rehab and rehabilitation therapy

HYGIENE



Shower and commode chairs facilitate a dignified, multi-positional hygiene experience

TRANSPORT



All the benefits of a standard Broda wheelchair, certified for vehicular travel


BEHAVIORAL



Stationary chairs with safe gliding that promote calm and contentment

Make Sure

USER BENEFITS*



*Users obtain these benefits when seated in a correctly-fitting Broda wheelchair with proper postural alignment.

THE BRODA DIFFERENCE

Comfort for Care



Comfort is essential, especially for those who spend much of their time in a wheelchair. Make comfortable seating part of long-term care with Broda's **ready out-of-the-box** wheelchairs.


Proven for Providers

Broda wheelchairs are meticulously designed to alleviate stress and injuries for caregivers and healthcare providers. With easy adjustments and customizable options, Broda products offer **best-in-class** functionality.



It's Broda.

CAREGIVER & FACILITY BENEFITS



Durable by Design



With over four decades of experience (est. 1981) and an **industry-leading**, limited lifetime frame warranty on 16-gauge steel, Broda provides safe, comfortable seating that stands the test of time.

Broda wheelchairs are trusted by patients, families, and caregivers to improve each patient's physical and emotional wellbeing, leading to **better quality of life**.

Time-Tested for Trust



SEATING & POSITIONING ESSENTIALS

Features and functionality for most Broda positioning, pedal, and rehab wheelchairs

SEAT TILT

Keeps the occupant's hips at a fixed angle

- Facilitates pressure redistribution
- Allows posterior positioning for comfort
- Angle-adjustable seat accommodates lifts and transfers



BACK RECLINE

Opens angle of the occupant's hips

- Accommodates hip flexion limitations
- Decreases postural deviations, such as lateral lean and slumping
- Trendelenburg-capable for maximum comfort

ADJUSTABLE & REMOVABLE ARM SUPPORTS

Provides support for the occupant's arms

- Improves postural stability with support of upper extremities
- Creates easy patient access for frequent repositioning
- Gives additional support for conditions such as edema



ADJUSTABLE & REMOVABLE LEG SUPPORTS


Creates support for the occupant's legs and/or feet

- Improves comfort and circulation lower extremity support
- Provides stability and improves posture
- Facilitates safe, obstacle-free, front standing transfers

ADJUSTABLE & REMOVABLE LATERAL SHOULDER SUPPORTS


Helps maintain the occupant's postural support

- Promotes a midline posture and minimizes postural deviations
- Easily removable for safe transfers
- Enhances users' comfort and safety




WHAT PEOPLE ARE SAYING ABOUT BRODA


Family Members

 At the age of 82, my father had advanced Parkinson's and dementia. His power wheel chair was taken away for safety concerns. In a fortunate twist of fate, he qualified for the VA to fund a Broda wheelchair. This wheelchair has improved his quality of life and I can't imagine him not being in it over the past few years. We are fortunate that the VA will fund Broda chairs for veterans like him who need them.


Family Member
Kansas, USA

 My son has had a Broda chair for the past 8 years. He was 23 years old when he acquired his [Traumatic Brain Injury] from a motorcycle accident. The Broda chair has allowed him to be up out of bed on a daily basis. If not for his Broda, he wouldn't be able to get out of bed because he has no truck control and cannot sit in other chairs for any extended periods.


Judy
Michigan, USA

 My brother who has ALS is completely paralyzed. We had a dream of taking him to his summer home in Vermont. The Synthesis [Positioning Wheelchair] went above and beyond to provide him with a comfortable wheelchair that allowed him to enjoy everything he wanted, including visiting with friends, hanging out at our boathouse, and even going on the boat itself! It was a nice change of pace over having to stay in bed all day.

Amy T.
California, USA




Clinicians

 Our staff really likes the Broda chairs. They are easy to use, adjust, roll, etc. They are excellent, durable chairs and we are very satisfied.

Kim K.
Georgia, USA

 My facility had been trying for months to accommodate our Huntington's resident with comfortable seating...As soon as he was placed in the Broda Elite Positioning Wheelchair, he never made a sound of pain or cried out in anger, and he actually tried to communicate with those around him. We were so pleased this gentleman did not need as much supervision and care as he had in the other chairs he had tried. We were so amazed at how quickly his behavior improved with Broda.

Physical Therapist
New York, USA

 Our facility had a gentleman with severe contractures in his legs and he had not been out of bed for a very long time...After discussing his needs, we arranged a trial of an Elite Positioning Wheelchair for this resident. For the first time in a very long time, the resident was able to be out of bed, comfortable, and social again. We want to thank Broda for solving an issue that we thought was unsolvable.

Occupational Therapist
South Carolina, USA



BRODA BENEFITS

At Broda, we promise to promote and preserve the wellbeing of the people who use our products. You can be sure your Broda chair or wheelchair will provide you with the following benefits.



SEATING COMFORT

We believe comfort is essential to a high quality of life and overall wellbeing.

Our products offer postural alignment, pressure relief, and secure seating.



POSITIONING



MOBILITY

Broda will help you reach the independence and social connection you deserve.

We'll build a wheelchair to match your style and provide the perfect fit.



PERSONALIZATION



Each of our products protects the user with built-in safety features.

PATIENT SAFETY

We are dedicated to making healthcare facilities safe for every caregiver.



CAREGIVER SAFETY



Keeping you healthy is our priority. Our products are easy to clean and disinfect.

INFECTION CONTROL

Our long-lasting products are guaranteed by our unbeatable 10-year frame warranty.




DURABILITY



Our team of seating and positioning experts are here to guide you through the process.

EASY ORDERING

POSITIONING

 The **Synthesis Positioning Wheelchair** has been a blessing for my father, his caregivers and my family. It **improves his quality of life** and I can't imagine him not being in it over the past few years.

- Kent P.
Family Member



Synthesis *Positioning Wheelchair*

A FULL RANGE OF POSITIONING FEATURES
WITH CONTEMPORARY DESIGN



Shown with rear 14" mag wheel option

PRODUCT FEATURES

DEGREES OF TILT	Up to 40°
DEGREES OF RECLINE	Up to 90° & Trendelenburg Capable
SEAT WIDTHS	16", 18", 20", 22", 24", 26"
SEAT DEPTHS	17", 18.5", 20"
WEIGHT CAPACITY	Up to 350 lbs.



Tilt functionality provides additional pressure offloading



Front pivot seat tilt maintains a forward line of sight, enhancing engagement



Adjustable lateral shoulder supports fit a variety of users



Easy side transfers for safe patient handling

Elite *Positioning Wheelchair & Tilt Wheelchair*

ENHANCED POSITIONING AND FUNCTIONALITY



Shown with rear-mounted mag wheel option

PRODUCT FEATURES

ONLY AVAILABLE IN CANADA

DEGREES OF TILT
DEGREES OF RECLINE
SEAT WIDTHS
SEAT DEPTHS
WEIGHT CAPACITY

Positioning: Up to 40° | **Tilt:** Up to 40°
Positioning: Up to 90° | **Tilt:** N/A
16", 18", 20", 22", 24", 26", 28"
15.5", 17", 18.5", 20"
Up to 350 lbs.



Tilt and recline functions enhance pressure redistribution and postural support



*Accommodates **safe patient handling** and changing of incontinence products*



Removable arm supports make transfers safer for caregivers and patients



WC19 Transport Package available for safe patient vehicle transport

Midline *Positioning Wheelchair*

*ENHANCED POSITIONING
WITH REFINED FEATURES*



Shown with rear-mounted mag wheel option

PRODUCT FEATURES

DEGREES OF TILT	Up to 38°
DEGREES OF RECLINE	Up to 90°
SEAT WIDTHS	16", 18", 20", 22", 24", 26", 28"
SEAT DEPTHS	15.5", 17", 18.5", 20"
WEIGHT CAPACITY	Up to 350 lbs.




Full tilt and recline provides superior repositioning



Adjustable seat surface height adapts the Midline to support users at multiple care levels



Tilt functionality enhances comfort



Comfort Tension Seating® increases support as well as sitting tolerance

Vanguard *Bariatric Wheelchair*

*BARIATRIC POSITIONING
CARE REDEFINED*



Shown with the Acute Care Configuration (ACC) full vinyl wrap option


PRODUCT FEATURES

*DEGREES OF TILT
DEGREES OF RECLINE
SEAT WIDTHS
SEAT DEPTHS
WEIGHT CAPACITY*

*Up to 30° Posterior / Up to 10° Anterior
Up to 45°
24", 26", 28", 30"
Adjustable 21" or 23"
Up to 600 lbs.*




Anterior and posterior tilt aids in safe patient handling



Swing away, removable arm supports facilitate side-transfers



Manual foot pump allows safe, no-lift tilt positioning



Optional concave back for gluteal shelf improves user support & dignity

PEDAL



Prior to the **Encore wheelchair**, we had tried numerous positioning devices for the resident, which were ineffective. She had a tendency to lean forward and to one side. Once using the Encore, the resident's posture immediately improved. She, herself, was very thankful because she was "**so comfortable**". The staff couldn't believe the result. Besides this resident sitting in proper alignment, she could still self-propel with her feet.

- Linda P.
Restorative Coordinator

Encore *Pedal Wheelchair*

MOBILITY WITH COMFORT,
FUNCTIONALITY, & STYLE



HME B
New Product
Award 2019

TODAY'S
CAREGIVER
2019
FRIENDLY
AWARD
CAREGIVER.COM

mobility
product award 2019



Dynamic frame design includes
adjustable shoulder supports
and tall back option



Front pivot seat tilt optimizes
hand and foot propulsion



Optional dynamic rocking calms
and creates contentment

Shown with short back and adjustable head support

PRODUCT FEATURES

DEGREES OF TILT	Up to 40°
SEAT WIDTHS	16", 18", 20", 22", 24", 26", 28" 16", 18", 20", 22" with Dynamic Rocking
SEAT DEPTHS	17", 18.5", 20" 15.5", 17", 18.5", 20" with Dynamic Rocking
WEIGHT CAPACITY	350 lbs., 250 lbs. with Dynamic Rocking
SEAT HEIGHTS	12.5", 13.5", 14.5", 15.5", 16.5", 17.5" 14", 15", 16", 17", 18", 19" with Dynamic Rocking



Optional WC19 transport package
for secure vehicle travel

Sashay & Latitude

Pedal Wheelchairs

ENHANCED MOBILITY AND FREEDOM

SASHAY PEDAL WHEELCHAIR

LATITUDE PEDAL WHEELCHAIR



SASHAY FEATURES

LATITUDE FEATURES

DEGREES OF TILT
SEAT WIDTHS
SEAT DEPTHS
WEIGHT CAPACITY
SEAT HEIGHTS

Up to 18°
16", 18", 20", 22", 24", 26", 28"
15.5", 17", 18.5", 20"
Up to 350 lbs.
12.5", 13.5", 14.5", 15.5", 16.5", 17.5"

Up to 33° with Dynamic Rocking
16", 18", 20", 22"
15.5", 17", 18.5", 20"
Up to 250 lbs.
14", 15", 16", 17", 18", 19"



Sashay **low seat surface height** facilitates foot propulsion



Latitude **front mount mag wheels** enable hand propulsion



Latitude **gentle rocking motion** soothes agitation and creates contentment

REHABILITATION (E1161)



Broda's wheelchairs make it **easier for caregivers** to provide optimal care with less stress on the client and the caregiver. For example, easier transfers, easier re-positioning, less maintenance, and easier mobility.

- Carole P.
MS/OTR.

Encore *Rehab Wheelchair*

MOBILITY WITH COMFORT,
FUNCTIONALITY, &
CUSTOMIZED SEATING



Shown with Matrix Libra cushion, Matrix E2 back support, Matrix Elan head support



Back canes fold down for easy storage options



Dynamic rocking function soothes and stimulates users



Quick-release mag wheels can be removed for easy transport



Optional WC19 Transport Package for secure vehicle travel

PRODUCT FEATURES

DEGREES OF TILT

Up to 33° with Dynamic Rocking

SEAT WIDTHS

16", 18", 20", 22"

SEAT DEPTHS

16", 17", 18", 19", 20"

SEAT HEIGHTS

Adjustable 14", 15", 16", 17", 18", 19"

WEIGHT CAPACITY

Up to 250 lbs.

Synthesis *Rehab Wheelchair*

*DYNAMIC POSITIONING
THROUGHOUT REHABILITATION*



HCPCS
E1161




Tilt functionality offloads pressure and promotes upright posture



Front pivot seat tilt maintains a forward line of sight, enhancing engagement



Adjustable upper lateral supports fit a variety of users



Solid seat pan provides support for a variety of cushions

Shown with rear 20" mag wheels and Matrix PS Cushion

PRODUCT FEATURES

DEGREES OF TILT	Up to 40°
DEGREES OF RECLINE	Up to 90°
SEAT WIDTHS	16", 18", 20", 22", 24", 26"
SEAT DEPTHS	17", 18.5", 20"
SEAT HEIGHTS	19.5", 20.5", 22", 23.5"
WEIGHT CAPACITY	Up to 350 lbs.

Comfort *Rehab Wheelchair*

SIMPLIFIED DURABILITY FOR REHABILITATION CARE



HCPCS
E1161

Shown with Matrix PS cushion, Matrix Elite TR back support, Matrix Elan Head Support, and swing-away lower leg supports

PRODUCT FEATURES

DEGREES OF TILT	Up to 45°
SEAT WIDTHS	16", 18", 20", 22"
SEAT DEPTHS	Adjustable 16", 17", 18", 19", 20"
SEAT HEIGHTS	Adjustable 14", 15", 15.5", 16", 16.5", 17", 17.5", 18", 19"
WEIGHT CAPACITY	Up to 250 lbs.



Front pivot seat tilt maintains user comfort with smooth, effortless tilt



Posterior tilt optimizes pressure redistribution



Solid seat pan offers support and is compatible with a variety of seat cushion options



Range of simple adjustments accommodates a variety of needs

Elite Rehab Wheelchair

ENHANCED MOBILITY
AND FREEDOM



Shown with rear 20" mag wheels and Matrix PSP cushion

HCPCS
E1161




Solid seat pan option
accommodates a variety of Matrix
cushions for additional support



Adjustable upper lateral supports
fit a variety of users



Posterior tilt effectively opens the
diaphragm improving digestion,
oxygenation, and organ function



Optional Comfort Tension Seating®
improves postural support and
sitting tolerance

PRODUCT FEATURES

DEGREES OF TILT	Up to 40°
DEGREES OF RECLINE	Up to 15°
SEAT WIDTHS	16", 18", 20", 22", 24", 26", 28"
SEAT DEPTHS	17", 18.5", 20"
SEAT HEIGHTS	20.5"
WEIGHT CAPACITY	Up to 350 lbs.

TRANSPORT



For facilities, we are **using Broda to reduce their cost of patient transportation**. For the consumers it's just comfort. Broda provides 100% comfort; you can't even compare the difference between Broda wheelchairs to other wheelchairs.



- Weslava P.
Graceful Living Medical Transportation

Synthesis *Transport Wheelchair*

COMFORT AND SAFETY
ON THE GO



Shown with rear 14" mag wheels

PRODUCT FEATURES


DEGREES OF TILT	Up to 40°
DEGREES OF RECLINE	Up to 90°
SEAT WIDTHS	16", 18", 20", 22", 24"
SEAT DEPTHS	17", 18.5", 20"
WEIGHT CAPACITY	Up to 350 lbs.




Crash tested and WC-19 compliant for vehicle transport



Easy to load in a variety of vehicles



Brackets enable safe patient transport between locations



Tilting functionality allows for adjustable comfort in transit

HYGIENE



Broda's **Revive Shower Commode** meets our patient's needs. He indicates he feels safe and he is able to participate in his own care each morning. He is **safe at all times** while in the chair, and the caregivers are safe at all times while providing his care. All staff report being pleased with the commode and find it much easier to assist him.

- Pam
Director of Nursing

Revive Shower Commode

*DURABLE BATHING SYSTEM
FOR SAFE PATIENT CARE*



Anterior and posterior tilt aid in fall prevention and safe patient handling



Multi-positional functionality with tool-free adjustability for ease of use



Multiple configurations available to accommodate unique needs



Rust-resistant frame and components, with a 10-year warranty

PRODUCT FEATURES

<i>DEGREES OF TILT</i>	<i>Up to 12° Posterior, Up to 10° Anterior</i>
<i>DEGREES OF RECLINE</i>	<i>Up to 67°</i>
<i>SEAT WIDTH</i>	<i>18", 22"</i>
<i>SEAT DEPTH</i>	<i>18"</i>
<i>WEIGHT CAPACITY</i>	<i>Up to 400 lbs.</i>

Sentinel Shower Commode

SAFE, COMFORTABLE,
INDEPENDENT BATHING



Shown with standard mag wheel configuration and optional commode pan and basket

PRODUCT FEATURES

BACK ANGLE	Adjustable 5°, 10°, 15°
SEAT ANGLE	Adjustable 0° - 6°
SEAT WIDTH	18", 20", 22"
SEAT DEPTH	Adjustable 16 - 19"
WEIGHT CAPACITY	Up to 300 lbs.



Features a **range of simple adjustments** for seat depth, wheel position, seat height, and more.



2-in-1 solid seat and commode is highly-versatile and easily mountable in all four directions



Configurable with either **user-propelled mag wheels** or cost-effective casters



Flip-up arm supports and swing-away leg supports allow for easy bathing and transfers

Revive *Bariatric Shower Commode*

*BARIATRIC HYGIENE CARE IN
A CLASS ALL ITS OWN*



PRODUCT FEATURES

DEGREES OF TILT	<i>Up to 14° Posterior, Up to 8° Anterior</i>
DEGREES OF RECLINE	<i>Up to 67°</i>
SEAT WIDTH	<i>27", 33" with width adjustable armrests</i>
SEAT DEPTH	<i>21", 23"</i>
WEIGHT CAPACITY	<i>Up to 800 lbs.</i>



Anterior tilt aids in patient transfers and ease of patient care



Mid-mounted Mag Wheel provides increased maneuverability



Optional concave back accommodates and supports gluteal shelf



Rust-resistant frame and components are durable and long-lasting

BEHAVIORAL



*Our Tranquille gliders are always disappearing out of our common areas and end up in our patients' rooms. We don't know how they get there, but they love them and **prefer sitting in them** more than anything else. One patient with a back injury told us this was the only chair he could sit comfortably in. He wanted to take it home with him after his rehabilitation was complete.*

*- Josiah D.
Facility Administrator*

Tranquille & Aspire Gliders


TRANQUILLE GLIDER



ASPIRE GLIDER




Gentle **gliding motion**



Stimulates active patients and reduces the need to wander



Auto-locking feature increases safety and prevents falls



Comfort Tension Seating® for immersive, all-day comfort

TRANQUILLE

ASPIRE

SEAT WIDTHS

18", 20"

12", 14", 16"

SEAT DEPTHS

17"

10", 12", 14"

BACK LENGTHS

26", 32"

18", 20", 22"

WEIGHT CAPACITY

Up to 250 lbs.

Up to 200 lbs.

Padding Packages

ADDITIONAL POSITIONING PADDING (APP)



- Seat Cushion
- Back and Head Support
- Foot and Calf Support
- Lateral Shoulder Support Cushions

APP shown on the Synthesis Positioning Wheelchair

RECOMMENDED FOR:

Use with positioning chairs to give patients extra comfort.

SLIM POSTURE PADDING (SPP)



- Seat Cushion
- Back Support
- Extended Back Support*
*for tall back height options

SPP shown on the Encore Pedal Wheelchair

RECOMMENDED FOR:

Use with pedal chairs for extra comfort with a lower profile.

HUNTINGTON'S SPECIAL PADDING (HSP)



- Long Lateral Shoulder Supports
- Thigh Positioning Strap
- Long Footbox

HSP shown on the Synthesis Positioning Wheelchair

RECOMMENDED FOR:

Users with neurological and movement disorders.

ACUTE CARE CONFIGURATION (ACC)



- Full Vinyl Covers
- ABS Footrest with Calf and Sole Pads
- Covers Comfort Tension Seating® straps for easy cleaning

ACC shown on the Synthesis Positioning Wheelchair

RECOMMENDED FOR:

Healthcare environments where infection control and durability are a high priority.

Rehab Seating

matrx[®] SEATING SERIES

For convenience, all Matrx[®] options ship on-chair for Broda Rehab Wheelchairs

LIBRA CUSHION



VI CUSHION



PSP CUSHION



PS CUSHION



ELITE E2™ BACK



ELITE BACKS (TR/DEEP)



PB BACK (DEEP)



*LOXX
HEAD SUPPORT*



*ELAN
HEAD SUPPORT*

The Matrx[®] Seating series pairs perfectly with Broda's Rehab wheelchairs. Additional cushion, back, and headrest accessories available.

Accessories

Call us for additional accessory options.
1-844-552-7632

FRAME ACCESSORIES

- O² Soft Carry Bag
- IV Pole
- Rear, Front, or Lateral Stabilizers
- Lateral Torso Supports



O² Soft Carry Bag



Lateral Supports

PADDING AND ARM SUPPORTS

- Long Arm Supports
- Side Arm Support
- Side Reduction Pads
- Arm Trough



Arm Trough



Side Arm Support

TRAYS

- Slide-On Full Tray
- Slide-On Half Tray
- Padded Tray Cover



Slide-On Half Tray



Slide-On Full Tray

LEG & FOOT SUPPORTS

- ABS Leg Support
- Leg Support Panel
- Foot Box
- Leg Support
- Swing-Away Leg Support



ABS Legrest



Foot Box

HEAD & NECK SUPPORTS

- Contoured Neckrest Pad
- Neckrest with Two Round Bolsters



Neck Support with Bolsters



Contoured Head Support

POSITIONING STRAPS

- Thigh Positioning Belt
- Pelvic Positioning Belt



Padded Thigh Belt



Pelvic Positioning Belt

Chair Specifications

Available Only in Canada



Synthesis Positioning Wheelchair
Synthesis Transport Wheelchair



Elite Positioning Wheelchair
Elite Tilt Wheelchair



Midline Positioning
Wheelchair



Centric Positioning
Wheelchair



Vanguard® Bariatric Wheelchair

Weight Capacity	Up to 350 lbs.	Up to 350 lbs.	Up to 350 lbs.	Up to 250 lbs.	Up to 600 lbs.
Degrees of Tilt	Up to 40°	Up to 40°	Up to 38°	Up to 25°	Up to 30° Posterior Up to 10° Anterior
Degrees of Recline	Up to 90° - Trendelenburg capable	Positioning: Up to 90° Tilt: N/A	Up to 90° - Trendelenburg capable	Up to 45°	Up to 45°
Seat Widths Available	16", 18", 20", 22", 24", 26"	16", 18", 20", 22", 24", 26", 28"	16", 18", 20", 22", 24", 26", 28"	16", 18", 20", 22", 24", 26", 28"	24", 26", 28", 30"
Seat Depths	Positioning: 17", 18.5", 20", 22" Transport: 17", 18.5", 20"	15.5", 17", 18.5", 20", 22"	15.5", 17", 18.5", 20"	15.5", 17", 18.5", 20"	Adjustable 21" or 23"
Back Heights	32", 33.5", 35"	32", 33.5", 35"	32", 33.5", 35"	32", 33.5", 35"	Adjustable 32", 34", 36", 38", 40"
Seat Surface Heights	19.5", 20.5", 22", 23.5"	19.5", 20.5", 22", 23.5"	Adjustable 17", 18", 19", 20", 21", 22"	Adjustable 19.5", 20.5", 21.5"	Adjustable 20", 21", 22", 23", 24"
Mag Wheels	14" Rear Mounted 14" Mid Mounted 20" Rear Mounted	14" Mid Mounted 14" Rear Mounted 20" Rear Mounted	20", 22", 24" Front Mounted 20", 22" Rear Mounted	20", 22", 24" Front Mounted 20", 22", 24" Rear Mounted	N/A
Recommended Seating	Comfort Tension Seating® Additional Positioning Padding (APP) Cushion	Comfort Tension Seating® Additional Positioning Padding (APP) Cushion	Comfort Tension Seating® Additional Positioning Padding (APP) Cushion	Comfort Tension Seating® Additional Positioning Padding (APP) Cushion	Comfort Tension Seating®
Recommended Back Support	Comfort Tension Seating® Additional Positioning Padding (APP) Back	Comfort Tension Seating® Additional Positioning Padding (APP) Back	Comfort Tension Seating® Additional Positioning Padding (APP) Back	Comfort Tension Seating® Additional Positioning Padding (APP) Back	Comfort Tension Seating®
Recommended Head Support	Broda® Neck Support	Broda® Neck Support	Broda® Neck Support	Broda® Neck Support	Broda® Neck Support

Chair Specifications



Encore Pedal Wheelchair



Sashay Pedal Wheelchair



Latitude Pedal Wheelchair



Revive Shower Commode
Revive Bariatric
Shower Commode



Sentinel
Shower Commode



Tranquille Glider
Aspire Pediatric Glider

Weight Capacity	Up to 350 lbs. Up to 250 lbs. (with rocking)	Up to 350 lbs.	Up to 250 lbs.	Up to 400 lbs. Bari: up to 800 lbs.	Up to 300 lbs.	Up to 250 lb. Aspire: up to 200 lb.
Degrees of Tilt	Up to 40°	Up to 18°	Up to 33°	Up to 12° Posterior Up to 10° Anterior Bari: Up to 14° Posterior Up to 8° Anterior	5°, 10°, 15°	N/A
Degrees of Recline	N/A	N/A	N/A	Up to 67°	N/A	N/A
Seat Widths	16", 18", 20", 22", 24", 26" 28" 16", 18", 20", 22" (with dynamic rocking)	16", 18", 20", 22", 24", 26", 28"	16", 18", 20", 22"	18", 22" Bari: 27"-33"	18", 20", 22"	18", 20" Aspire: 12", 14", 16"
Seat Depths	17", 18.5", 20"	15.5", 17", 18.5" or 20"	15.5", 17", 18.5" or 20"	18" Bari: 21", 23"	Adjustable 16" - 19"	17" Aspire: 10", 12", 14"
Back Heights	22", 32", 33.5", 35"	32", 33.5", 35"	32", 33.5", 35"	26"	Adjustable 16", 18", 20"	26", 32" Aspire: 18", 20", 22"
Seat Surface Heights	Adjustable 12.5" - 17.5" Adjustable 14" - 19" (with dynamic rocking)	Adjustable 12.5", 13.5", 14.5", 15.5, 16.5", 17.5"	Adjustable 14", 15", 16", 17, 18" 19"	24.5"	19.5, 23.5	16", 18"
Mag Wheels	20" or 22" Front Mounted	20", 22" Front Mounted 20", 22" Rear Mounted	20", 22" Front Mounted 20", 22" Rear Mounted	N/A	24" Rear Mounted	N/A
Recommended Seating	Comfort Tension Seating* Slim Positioning Padding (SPP) Cushion	Comfort Tension Seating* Slim Positioning Padding (SPP) Cushion	Comfort Tension Seating* Slim Positioning Padding (SPP) Cushion	N/A	N/A	N/A
Recommended Back Support	Comfort Tension Seating* Slim Posture Padding (SPP) Back	Comfort Tension Seating* Slim Posture Padding (SPP) Back	Comfort Tension Seating* Slim Posture Padding (SPP) Back	N/A	N/A	N/A.
Recommended Head Support	Matrx® Elan Headrest Broda® Neck Support (tall back)	Broda® Neck Support	Broda® Neck Support	Broda® Neck Support Broda® Head Support	Broda® Neck Support Broda® Head Support	Broda® Neck Support

Chair Specifications



Encore Rehab Wheelchair



Comfort Rehab Wheelchair



Synthesis Rehab Wheelchair



Elite Rehab Wheelchair



Latitude Rehab Wheelchair

Weight Capacity	Up to 250 lbs.	Up to 250 lbs.	Up to 350 lbs.	Up to 350 lbs.	Up to 250 lbs.
Degrees of Tilt	Up to 33° (with dynamic rocking)	Up to 45°	Up to 40°	Up to 40°	Up to 33° (with dynamic rocking)
Degrees of Recline	N/A	N/A	Up to 90°	15°	N/A
Seat Widths	16", 18", 20", 22"	16", 18", 20", 22"	16", 18", 20", 22", 24", 26"	16", 18", 20", 22", 24", 26", 28"	16", 18", 20", 22"
Seat Depths	16", 17", 18", 19", 20"	16", 17", 18", 19", 20"	17", 18.5", 20"	17", 18.5", 20", 22"	15.5", 17", 18.5", 20"
Back Cane Angles	5°, 10°, 15°, 20°	5°, 10°, 15°, 20°	N/A	N/A	N/A
Back Heights	16", 18", 21", 24"	16", 18", 21", 24"	32", 33.5", 35"	32"	32", 33.5", 35"
Seat Surface Heights	Adjustable 14", 15", 16", 17", 18", 19"	Adjustable 14", 15", 15.5", 16", 16.5", 17", 17.5", 18", 19"	19.5", 20.5", 22", 23.5"	20.5"	Adjustable 14", 15", 16", 17", 18", 19"
Mag Wheels	20", 22" Front Mounted	20", 22", 24" Rear Mounted	14", 20" Rear Mounted 14" Mid Mounted	14", 20", 22" Rear Mounted	20", 22" Front or Rear Mounted
Recommended Seating	Matrx® Libra Cushion	Matrx® Libra Cushion	Matrx® Libra Cushion	Matrx® Libra Cushion	Matrx® Libra Cushion
Recommended Back Support	Matrx® E2 Back Support	Matrx® E2 Back Support	Matrx® E2 Back Support	Comfort Tension Seating® Slim Posture Padding (SPP) Back	Comfort Tension Seating® Slim Posture Padding (SPP) Back
Recommended Head Support	Matrx® Loxx Headrest	Matrx® Loxx Headrest	Broda Neck Support	Broda Neck Support	Broda Neck Support

Frame Options

BLACK	BRONZE	HAMMER GREY	TUNGSTEN SAND	SATIN NICKEL	CARMINE RED	OCEAN BLUE	ZINC YELLOW	DEEP ORANGE	PASSION PINK	ROYAL PURPLE
-------	--------	-------------	---------------	--------------	-------------	------------	-------------	-------------	--------------	--------------

	BLACK	BRONZE	HAMMER GREY	TUNGSTEN SAND	SATIN NICKEL	CARMINE RED	OCEAN BLUE	ZINC YELLOW	DEEP ORANGE	PASSION PINK	ROYAL PURPLE
Synthesis Positioning Wheelchair	✓	✓		✓	✓	✓	✓	✓		✓	✓
Synthesis Transport Wheelchair	✓	✓		✓	✓	✓	✓	✓		✓	✓
Midline Positioning Wheelchair	✓		✓								
Elite Tilt Wheelchair	✓		✓								
Elite Positioning Wheelchair (CAN)	✓		✓								
Centric Positioning Wheelchair			✓								
Vanguard Positioning Wheelchair			✓								
Encore Pedal Wheelchair	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓
Latitude Pedal Wheelchair	✓		✓								
Sashay Pedal Wheelchair	✓		✓								
Encore Rehab Wheelchair	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓
Comfort Rehab Wheelchair	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓
Synthesis Rehab Wheelchair	✓	✓		✓	✓	✓	✓	✓		✓	✓
Elite Rehab Wheelchair	✓		✓								
Latitude Rehab Wheelchair	✓		✓								
Tranquille Glider		✓	✓	✓	✓	✓	✓				
Aspire Glider		✓		✓	✓	✓	✓				

Upholstery Options



	SAVOY BLACK	INTERLACE GRAPHITE	INTERLACE ESPRESSO	INTERLACE INDIGO	BEELINE FORGE	BEELINE PORTOBELLO
Synthesis Positioning Wheelchair	✓	✓	✓	✓	✓	✓
Synthesis Transport Wheelchair	✓	✓	✓	✓	✓	✓
Midline Positioning Wheelchair						
Elite Tilt Wheelchair						
Elite Positioning Wheelchair (CAN)						
Centric Positioning Wheelchair						
Vanguard Positioning Wheelchair						
Encore Pedal Wheelchair	✓	✓	✓	✓	✓	✓
Latitude Pedal Wheelchair	✓	✓	✓	✓	✓	✓
Sashay Pedal Wheelchair	✓	✓	✓	✓	✓	✓
Encore Rehab Wheelchair						
Comfort Rehab Wheelchair						
Synthesis Rehab Wheelchair	✓	✓	✓	✓	✓	✓
Elite Rehab Wheelchair						
Latitude Rehab Wheelchair	✓	✓	✓	✓	✓	✓
Tranquille Glider	✓	✓	✓	✓		✓
Aspire Glider				✓		

Upholstery Options



	BEELINE GEYSER	BEELINE FOSSIL	SILVERTEX CARBON	SILVERTEX MOCHA	SILVERTEX SAPPHIRE	SILVERTEX ACID
Synthesis Positioning Wheelchair	✓	✓				
Synthesis Transport Wheelchair	✓	✓				
Midline Positioning Wheelchair					✓	
Elite Tilt Wheelchair					✓	
Elite Positioning Wheelchair (CAN)					✓	
Centric Positioning Wheelchair					✓	
Vanguard Positioning Wheelchair					✓	
Encore Pedal Wheelchair	✓					
Latitude Pedal Wheelchair	✓	✓			✓	
Sashay Pedal Wheelchair	✓	✓			✓	
Encore Rehab Wheelchair						
Comfort Rehab Wheelchair						
Synthesis Rehab Wheelchair	✓	✓				
Elite Rehab Wheelchair					✓	
Latitude Rehab Wheelchair	✓	✓			✓	
Tranquille Glider	✓		✓	✓	✓	✓
Aspire Glider	✓	✓	✓	✓		✓


Upholstery Options



	RENAISSANCE BURGUNDY	RENAISSANCE METEOR	RENAISSANCE CHARCOAL	KALEIDOSCOPE CHOCOLATE	KALEIDOSCOPE SMOKEHOUSE	RETRO CAMO	DIGITAL CAMO
Aspire Glider		✓					
Centric Positioning Wheelchair	✓		✓				
Comfort Rehab Wheelchair							
Elite Positioning Wheelchair (CAN)	✓	✓	✓				
Elite Rehab Wheelchair	✓	✓	✓				
Elite Tilt Wheelchair	✓	✓	✓				
Encore Pedal Wheelchair					✓	✓	✓
Encore Rehab Wheelchair							
Latitude Pedal Wheelchair	✓		✓	✓	✓	✓	✓
Latitude Rehab Wheelchair	✓		✓	✓	✓	✓	✓
Midline Positioning Wheelchair	✓		✓				
Sashay Pedal Wheelchair	✓		✓	✓	✓	✓	✓
Synthesis Positioning Wheelchair				✓	✓	✓	✓
Synthesis Rehab Wheelchair				✓	✓	✓	✓
Synthesis Transport Wheelchair				✓	✓	✓	✓
Tranquille Glider	✓	✓	✓	✓	✓		
Vanguard Positioning Wheelchair	✓		✓				

Configure Your Wheelchair

Step 1: Select Chair Type



Step 2a: Select Positioning Wheelchair (Caregiver-Propelled)

Features, Functionality, and Customization Level

Bariatric

Vanguard Bariatric Wheelchair



Best Basic

 <p>Synthesis Positioning Wheelchair Synthesis Transport Wheelchair</p> <p>350 lbs. weight capacity Up to 40° Tilt Up to 90° Recline</p>	 <p>Elite Positioning Wheelchair Elite Tilt Wheelchair</p> <p>350 lbs. weight capacity Up to 40° Tilt Up to 90° Recline</p>	 <p>Midline Positioning Wheelchair</p> <p>350 lbs. weight capacity Up to 38° Tilt Up to 90° Recline</p>	 <p>Centric Positioning Wheelchair</p> <p>250 lbs. weight capacity Up to 25° Tilt Up to 45° Recline</p>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Optional Add-ons:

- Additional Positioning Padding (APP)
- Huntington's Special Padding (HSP)
- WC19 Transport Package
- Acute Care Configuration (ACC)
- Slim Posture Padding (SPP)

Step 2b: Select a Pedal Wheelchair - Self-Propelled (Hand or Foot)

Features, Functionality, and Customization Level

Best Basic

 <p>Encore Pedal Wheelchair with Dynamic Rocking feature</p> <p>250 lbs. weight capacity Up to 33° Tilt</p>	 <p>Encore Pedal Wheelchair</p> <p>350 lbs. weight capacity Up to 25° Tilt</p>	 <p>Latitude Pedal Wheelchair with Dynamic Rocking Feature</p> <p>250 lbs. weight capacity Up to 33° Tilt</p>	 <p>Sashay Pedal Wheelchair</p> <p>350 lbs. weight capacity Up to 18° Tilt</p>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Optional Add-ons:

- Slim Posture Padding (SPP)
- Additional Positioning Padding (APP)
- WC19 Transport Package

Step 3: Select Accessories, Colors, & Customizations

[See page 30 for common accessories.](#)

[See page 34 for frame and upholstery color options.](#)

Call us for ordering options.
1-844-552-7632

Configure Your Rehab Chair

More Independent

Less Independent

Step 2c: Select Chair Frame

BRODA



Encore Rehab Wheelchair

Latitude Rehab Wheelchair

Comfort Rehab Wheelchair

Synthesis Rehab Wheelchair

Elite Rehab Wheelchair

250 lbs. weight capacity Up to 33° Tilt Dynamic Rocking HCPCS E1161	250 lbs. weight capacity Up to 33° Tilt Dynamic Rocking HCPCS E1161	250 lbs. weight capacity Up to 45° Tilt Highly Adjustable HCPCS E1161	350 lbs. weight capacity Up to 40° Tilt Up to 90° Recline HCPCS E1161	350 lbs. weight capacity Up to 40° Tilt Up to 15° Recline HCPCS E1161
------------------------------------------------------------------------------	------------------------------------------------------------------------------	--------------------------------------------------------------------------------	--------------------------------------------------------------------------------	--------------------------------------------------------------------------------

For convenience, all Matrx® options ship on chair.

Step 3: Select Seating

matrx
SEATING SERIES

More Pressure Reduction

Less Pressure Reduction



Libra Cushion



Vi Cushion



PSP Cushion



PS Cushion

Encore Rehab or Comfort Rehab

Latitude Rehab or Elite Rehab

Step 4: Select Back Support

matrx
SEATING SERIES

More Supportive

Less Supportive



E2™ Backs



Elite (TR/Deep) Backs



PB & PB Deep Backs

HCPCS E2620/E2621	HCPCS E2620/E2621	HCPCS E2615	HCPCS E2601/E2611
-------------------	-------------------	-------------	-------------------

BRODA



*Comfort Tension Seating®
Additional Positioning Padding (APP)
Slim Posture Padding (SPP)*

Step 5: Select Head Support

matrx
SEATING SERIES

More Supportive

Less Supportive



Loxx head support
HCPCS E0955/E1028



Elan head support
HCPCS E0955/E1028

Step 6: Select Additional Accessories

Call us for more options!
1-844-552-7632

HCPCS codes supplied are for reference only and do not ensure coverage or payment. The HCPCS provided are based on our interpretation of the code definitions. HCPCS not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the M program and should consult an attorney or other advisor to discuss specific situations in further detail.



Broda USA
15455 Conway Rd, Suite 200
Chesterfield, MO 63017, USA

GSA Contract # 36F79721D0039

BRODA[®]
Essential Comfort

brodaseating.com
1-844-552-7632

Broda Canada
560 Bingemans Centre Dr
Kitchener, ON N2B 3X9, Canada

 **Ontario**
Assistive Devices Program