

BRODA[®]

Synthesis

Transport Wheelchair



Crash-tested and WC-19 compliant for vehicle transport



Easy to load in a variety of vehicles



Brackets enable safe patient transport between locations

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Promotes safe patient handling with enhanced positioning
- Proper seating alignment reduces patient fatigue and discomfort




Heavy-Duty Wheel Upgrade includes wheel locks and wide front casters

brodaseating.com | 1-844-552-7632

COMFORT AND SAFETY
ON THE GO


UPHOLSTERY OPTIONS



SPECIFICATIONS

Degrees of Tilt	Up to 40°
Degrees of Recline	Up to 90°
Weight Capacity	Up to 350 lbs. (158 kg)
Seat Widths	16", 18", 20", 22", 24" (40cm, 45cm, 50cm, 55cm, 63cm)
Seat Depths	17", 18.5", 20" (43cm, 46cm, 50cm)
Floor-to-Seat Heights	20.5" (52cm)
Wheels	5" (12cm) Aluminum Non-locking Front Casters 5" (12cm) Locking Twin Wheel Front Casters 6" (15cm) Aluminum Non-locking Front Casters 14" (35cm) Mid Mags 14" or 20" (35cm or 50cm) Rear Mags 6" (15cm) Multi-Directional Locking Rear Casters
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®
Standard Padding	Headrest, Shoulder Bolsters, Lower Side Pads, Seat & Back Pad with Cover

FRAME OPTIONS



Additional accessories available: foot supports, specialty padding packages, trays, and positioning straps